

Fully Integrated Spirometry released for AEROS

The fully integrated spirometry solution from InPS Alliance Partner Numed Healthcare, now allows Vision AEROS and LAN users to file high quality spirometry reports and Read coded data with ease. It is a very quick, simple and error free process that could save your nursing team around 10 minutes per spirometry test.

The TrueFlow spirometry solution utilises the latest ultrasonic flow measurement technology for supreme accuracy and durability. Numed provide a complete service that includes installation, product training and MediServe support to ensure that your service runs smoothly from the outset.



Key points

- Automatically files 15 spirometry Read codes and a full spirometry test report.
- Saves approximately 10 minutes administration time per spirometry test.
- Now available for LAN and AEROS Vision systems.

Advantages of full integration

1. The process is considerably faster than manual data entry and the reduced administration burden allows more patients to be tested in any given time.
2. Spirometry reports are filed directly to the patient's record as high quality PDFs. This makes them ideal for referral purposes and fully compatible with Choose and Book (to be replaced by NHS e-Referral in 2014) and GP2GP.
3. Spirometry Read codes and reports are always filed to the correct patient and there is a considerably reduced risk of data duplication using this integrated system.
4. All clinicians have immediate access to completed spirometry tests through Vision and hard copies can be printed at any Vision workstation within the practice if necessary.

What do existing Vision users think?

"Our Practice Nurse finds the new spirometer very easy to use - it also Read codes the numerous data automatically back onto the patient record saving her so much time and eliminating any chance of input error and likes the reports she can produce. The support and training have been quick and efficient"

Lisa Fogg, Practice Manager, Perranporth Surgery, Cornwall.


"We like the TrueFlow spirometer very much. It is easy to use with clear, easy to understand information on the blows appearing on screen. Excellent visual representation of an effective blow on screen with the balloon... a very useful aid for patients! Maintenance/cleaning is much easier than our old spirometer which saves time, as does the instant entry on to the Vision patient records. Definitely a good investment and any technical problems are very quickly and efficiently solved by your helpline technicians. Highly recommended!"

Rachel Young, Lead Practice Nurse, Sandown Health Centre, Isle of Wight.

How does it work?

The system is very simple to use:

1. The patient is selected within Vision, and their details are automatically transferred into the spirometry software.
2. The user can select which of the spirometry Read codes are to be filed to the patient's medical record at the end of the test, in addition to a PDF spirometry report:



The image shows the 'SpiroLink - Settings' window with the 'Read Codes' tab selected. It contains a section titled 'File record with the following Read Codes:' with checkboxes and input fields for various spirometry parameters. The parameters are grouped into columns: FVC, FEV1, PEF, SVC, and FEV1/FVC. Each parameter has a checkbox and a text input field. Below this, there are checkboxes for 'Add report comments to General Lung Function Test Read Code' and 'Patient Height (cm)', and a text input field for 'Patient Weight (kg)'. At the bottom, there is a button 'Enable Read Code Editing...' and 'Save' and 'Cancel' buttons.

	FVC	FEV1	PEF	SVC	FEV1/FVC
Actual	<input checked="" type="checkbox"/> 3396.	<input checked="" type="checkbox"/> 3390.	<input checked="" type="checkbox"/> 339A.	<input checked="" type="checkbox"/> 3386.	<input checked="" type="checkbox"/> 339R.
Predict	<input checked="" type="checkbox"/> 339Q.	<input checked="" type="checkbox"/> 339P.			
%	<input checked="" type="checkbox"/> 339i.	<input checked="" type="checkbox"/> 339S.			
Post	<input checked="" type="checkbox"/> 339h.	<input checked="" type="checkbox"/> 339b.	<input checked="" type="checkbox"/> 339B.		


Add report comments to General Lung Function Test Read Code ☒ 3371.

Patient Height (cm) ☒ 229.. Patient Weight (kg) ☒ 22A..

Enable Read Code Editing...

Save Cancel

3. The relaxed, forced and post bronchodilator (if required) spirometry tests are performed. Flow/Volume and Volume/Time graphs are displayed clearly on the computer screen in addition to the numerical results and a variety of other information including predicted values, test variance and interpretation.



4. Once the spirometry test is complete, it is filed complete with any user comments into the patient's record with a single button click. This automatically files all spirometry Read codes and a high quality PDF report (including graphs) to the patient record - perfect for electronic referral purposes.


14/03/14

gics & Adverse F
nd Reviews
ference
al History

ition Findings
sations
neous
t Results
stration Exam
on Clinic
entions
egisters

extension

us not recorded
otion
tion missing
ns Due in Nex...
/04/1991 o/d
anus/Polio ++ 2 2...












 Spirometry 10:34 - Automatic entry by SpiroLink software


[Show Less...](#)


He Lung function testing done -----

Parameter	L	Pre		Post	
		Pred	Best	%Pred	Best
FVC	4.38	5.23	119	5.12	
FEV1	3.83	4.27	111	4.10	
FEV1/FVC		0.84	0.82	97	0.80
PEF	L/min	478.80	479.20	100	537.34

* Indicates below LLN or significant post change.
 Predicted Set: ERS/ECCS
 System Interpretation Set: NICE(2010)
 Session Quality
 Pre D - Result not repeatable (FEV1 Var=0.79L (18.5%);
 FVC Var=0.55L (10.5%))
 Post D - Only one acceptable trial (FEV1 Var= -L (); FV
 C Var= -L ())
 System Interpretation
 Pre Normal Spirometry
 Post No interpretation, not enough acceptable maneuvers

 Peak flow rate after bronchodilation = 537.3 L/min
 Peak flow rate before bronchodilation = 479.2 L/min
 FEV1/FVC percent = 82 %
 FEV1 after bronchodilation = 4.1 L
 Percent predicted FEV1 = 111 %
 Forced expired volume in 1 second = 4.27 L
 Expected FEV1 = 3.83 L
 FVC after bronchodilation = 5.1 L
 FVC/Expected FVC percent = 119 %
 Forced vital capacity - FVC = 5.23 L
 Expected FVC = 4.38 L

 Spirometry 10:32 - Automatic entry by SpiroLink software



 Lung vital capacity = 5.04 L

[Show Less...](#)

He Lung function testing done Pre

Parameter	Pred	Best	%Pred
VC	4.41	5.04	114
VCex	4.41	5.04	114

Predicted Set: ERS/ECCS
 System Interpretation Set: NICE(2010)
 Session Quality
 Pre D - Only one acceptable trial (VC Var= -L ())

 Weight 80 kgs BMI: 24.6 O/E - weight
 Height 1.8 metres O/E - height

Sample Spirometry Report

Patient Information

Name
ID 123456
Age 43 (17/11/1967)
Height 180 cm
Weight 90 kg BMI 27.8
Gender Male
Ethnic Caucasian
Asthma Possible
Smoker No

Test Information

Test Type FVC (ex only)
Test Date 01/02/2011 17:09:40
Post Time 01/02/2011 17:11:05
Interpretation NICE(2010)
Predicted ERS/ECCS
Value Selection Best Value
BTPS (IN/EX) 1.12/1.02
User ID


Test Result

Your FEV1 / Predicted: 132%

Your Lung Age: 43

Parameter	Pred	LLN	Pre Best	Trial 1	Trial 2	%Pred	Post Best	Trial 2	Trial 1	%Pred	%Chg
VCmax [L]	-	-	5.73	5.73	5.12	-	6.02	6.02	5.88	-	5
FEV1 [L]	4.00	3.16	5.28	5.28	4.85	132	5.33	5.33	5.24	133	1
FVC [L]	4.91	3.91	5.73	5.73	5.12	117	6.02	6.02	5.88	123	5
FEV1/VCmax	-	-	0.921	0.921	0.947	-	0.884	0.884	0.892	-	-4
FEV1/FVC	0.794	0.676	0.921	0.921	0.947	116	0.884	0.884	0.892	111	-4
PEF [L/min]	561	442	756	756	730	135	723	700	723	129	-4
MMEF [L/s]	4.34	2.63	6.88	6.88	7.08	158	7.06	7.06	7.10	163	3
Session Quality	Pre		D - Result not reproducible (FEV1 Var=0.43L (8.2%); FVC Var=0.61L (10.7%))								
	Post		C (FEV1 Var=0.08L (1.6%); FVC Var=0.15L (2.4%))								
System Interpretation	Pre		Normal Spirometry								
	Post		Normal Spirometry								

Caution: Poor session quality. Interpret with care



For a product brochure or demonstration of the TrueFlow spirometer for AEROS or LAN Vision systems, please contact Numed Healthcare:

Call: **0114 2433896** Email: sales@numed.co.uk Web: www.numed.co.uk