

Packing Instructions

All equipment must be packaged correctly to ensure items are not damaged in transit.

The best way to achieve this is to use a Numed Transit Case.



You can purchase a Numed Transit Case by [clicking here](#)

Alternatively, please use a box making sure there is **adequate packaging** around the equipment. See below for examples.


Example 1



There should always be **sufficient packaging** around the equipment. Never place the patient leads on top of the machine, this can damage the display screen.

Add more packaging to divide items, for example by using either bubble wrap or packing paper.

Example 2



Padded envelopes should **never** be used to transport the machine itself, as they do not offer sufficient protection.

A transit case or box must be used. (see examples 1 and 2)



Example 3

Please ensure you have packaged your items correctly according to these guidelines before sending.

We **cannot accept liability** for items damaged in transit due to inadequate protection and/or being packaged in unsuitable casing.

Customers will be **liable for the cost** of repairing equipment damaged caused due to incorrect packaging being used.

numed
healthcare

If you have any questions,
then please contact us on **0114 243 3896**
or email: **info@numed.co.uk**