





# TrueFlow Spirometry Consumables

Fax Your Order: 0114 2433413 • Next Day Delivery on orders received before noon!

Tick to order    Quantity

	<p><b>Disposable Mouthpieces</b> A100973</p>	<p>Box of 200 <b>£38.85</b></p>	<input type="checkbox"/>	
	<p><b>Spirette M</b> <b>Multiple Patient Breathing Tube</b> <b>(Sterilisable) - A100829</b></p>	<p>Box of 50 <b>£50.40</b></p>	<input type="checkbox"/>	
	<p><b>Standard One Piece Spirette</b> <b>(Sterilisable) - A100974</b></p>	<p>Pack of 50 <b>£57.75</b></p>	<input type="checkbox"/>	
	<p><b>PeraSafe Sterilising Agent</b> <b>x5 16.2g Sachets - A100935</b></p>	<p>5 x 16.2g Sachets <b>£12.50</b></p>	<input type="checkbox"/>	
	<p><b>PeraSafe Sterilising Agent</b> <b>81g Container - A100594</b></p>	<p>81g Container <b>£12.50</b></p>	<input type="checkbox"/>	

For the TrueFlow Spirometer there is a choice of two mouthpiece options:

**Option 1: Spirette M (sterilisable) + Disposable Mouthpiece** - Following a spirometry test the Spirette M can be sterilised using a solution of PeraSafe. Each Spirette M can be sterilised and reused many times. The Disposable Mouthpiece portion is discarded after each test.

**Option 2: Standard One Piece Spirette (sterilisable)** - Following a spirometry test the Standard One Piece Spirette can be sterilised using a solution of PeraSafe. The Spirette can either be sterilised and reused many times or discarded after each test for the highest possible level of cross infection control.

Note:

Each new TrueFlow spirometer will be supplied with a starter pack of consumables consisting of:

- Disposable Mouthpieces – Box of 200
- Spirette M multiple patient breathing tube (sterilisable) - Box of 50
- PeraSafe Sterilising Agent (5 x 16.2g sachets OR 81g container)

Name: ..... Date: ..... Your Order Ref: .....

Delivery Address must include Postcode & Tel No:

Invoice Address if different:

.....  
.....

.....  
.....

Tel: .....

Signature: .....

Prices are exclusive of VAT & postage and were correct at time of printing.

Alliance House, Roman Ridge Road. Sheffield S9 1GB

**Tel: 0114 243 3896 Fax: 0114 243 3413 sales@numed.co.uk www.numed.co.uk**