

MicroMedical / Carefusion Spirometry Consumables 2010

Telephone Orders: 0114 2433896 • Fax Your Order: 0114 2433413

Next Day Delivery on orders received before noon!

Tick to order Quantity

	<p>3 Litre Calibration Syringe For use with all Numed Spirometers.</p>	<p>£295.00</p>	<input type="checkbox"/>	
	<p>Disposable Adult Mouthpieces with one-way valve</p> <ul style="list-style-type: none"> • High quality wax coated • Prevents the patient inhaling during a spirometry or peak flow test • Significantly reduces the risk of cross-infection between patients • Can be used with most makes of Spirometer / Peak Flow meter (standard 30mm diameter). 	<p>Box of 500 £65.00</p>	<input type="checkbox"/>	
	<p>Disposable Adult Mouthpieces</p> <ul style="list-style-type: none"> • High quality wax coated cardboard disposable mouthpieces • Can be used with most makes of Spirometer / Peak Flow meter (standard 30mm diameter). 	<p>Box of 500 £25.00</p>	<input type="checkbox"/>	
	<p>MicroGard Filter Maximum bacterial and viral filtration MicroGard® filter is a disposable barrier type Filter intended to protect both patient and equipment, by preventing the transmission of pathogens in droplets and aerosolised particles between the patient and the spirometer.</p>	<p>Pack of 25 £22.95 Pack of 100 £79.95</p>	<input type="checkbox"/> <input type="checkbox"/>	
	<p>Thermal Printer Paper 5 rolls per pack.</p>	<p>£22.50</p>	<input type="checkbox"/>	
	<p>PeraSafe Sterilising Solution Proven efficacy against viruses, bacteria, mycobacteria, fungi and spores. 81g makes up to 5 litres of solution.</p>	<p>£11.80</p>	<input type="checkbox"/>	

Name: Date: Your Order Ref:

Delivery Address must include Postcode & Tel No:

Invoice Address if different:

.....
.....

.....
.....

Tel:

Signature:

www.numed.co.uk

You can order by post: Numed Healthcare, Alliance House, Roman Ridge Road, Sheffield S9 1GB

Prices are exclusive of VAT & postage and were correct at time of printing.