

Micro 1

Handheld Spirometer



Based on an award winning design, Micro 1 is the latest generation of pocket spirometer. Designed specifically for situations where low cost precision spirometry measurements are required, the Micro 1 is extremely lightweight and portable. The large backlit graphical display and flat touch sensitive keypad ensures ease of use and functionality in this ergonomically designed spirometer.

Using the supplied PC configuration software, the Micro 1 spirometer can be configured to display the desired spirometry indices.



Features and Benefits

- Proven accuracy and reliability
- Choose up to 9 indices from FEV1, FEV6, FVC, FEV1/FEV6, FEV1/FVC, PEF, FEF25-75, FEF75 and FEF25
- Results can be displayed as % of predicted, Standard Deviation (SD) or Z Scores
- Quality prompts for each manoeuvre
- Quick patient test mode
- Backlit Monochrome LCD display
- Rechargeable battery (2 x AA NiMH) giving a battery life of 30 hours
- PC configuration software for user specific preferences
- Small, lightweight and portable
- Dimensions – 162 x 61 x 30mm
- Weight – 152gms

2 Year Warranty

MediServe Support

- 90 Days Included
- 1, 2 or 3 Year Contracts Available

Optional Extras

- 3 Litre Accuracy Check Syringe

Consumables

- One Way Safety Mouthpieces
- PeraSafe Sterilising Agent