




# MAC 600 ECG System

## A compact, lightweight and highly portable 12 Lead ECG system

MAC 600 is ideal for busy GP surgeries where a PC may not be readily available, or for community nurses requiring a compact and lightweight ECG system. It is supplied with all of the required software and accessories and a market leading 5 year warranty.



## Safer, Faster, Easier- I<sup>3</sup> Integration



MAC 600 is compatible with Numed's I<sup>3</sup> integration software\*, which automatically retrieves a patient's details from EMIS, TPP SystemOne or InPS Vision and files a PDF report and ECG Read code upon test completion to their medical record.

This offers a number of benefits to both you and your patients:



### Enhances Patient Care

Dramatically simplifies the ECG test procedure allowing you to concentrate on caring for your patients



### Improves Patient Safety

Ensures test results are always filed to the correct medical record, improving patient safety



### Improves Quality

Files a high quality ECG report rather than a poor quality scan, for a more accurate diagnosis



### Reduces Administration

Eliminates manual data entry, scanning and attachment of test results, giving you more time with your patients

See I<sup>3</sup> in action on our website: [www.numed.co.uk/i3](http://www.numed.co.uk/i3)

## What do our customers think?

*“The I<sup>3</sup> software enables this to be a fully integrated device placing results Read coded into the medical record. It saves time but also means that ECGs are filed to the correct patient record 100% of the time.”*

**Anwar Hussain - Practice Manager, Tower Hamlets**



## Peace of mind - MediServe Support

Faulty diagnostic equipment can be a major inconvenience, leading to delays in consultations and even the cancellation of appointments. MAC 600 is supplied with our unique MediServe support package\*\*, which includes:

- Unlimited technical support from a dedicated support team
- Full equipment repair cover including parts and labour
- Loan equipment when yours is being serviced or repaired
- Full annual service and calibration to manufacturer's specification

This will give you the peace of mind that your service will continue to run smoothly no matter what may happen.



## Confidence in your new equipment – Installation & Training

It is vital that you are confident in the use of your new diagnostic equipment. For a small fee our dedicated customer support team can:

- Install your new product at a date and time that is convenient for you and in as many consulting rooms as you require
- Ensure you get the most out of your new device through a wide range of training options, including training videos in our dedicated support section, remote webinar sessions and on-site training\*\*\* if required

## MAC 600 Additional Features:

- **GE Marquette 12SL**  
High quality interpretation as used in the majority of UK hospitals
- **Long Lasting Battery**  
Can perform up to 250 ECG tests on one charge
- **Hook Up Advisor**  
Allows you to resolve quality issues prior to ECG acquisition
- **5 Year Warranty**  
For years of trouble free operation

## What's inside the package?

- MAC 600 12 Lead ECG System
- 100 Ambu SUPAtab ECG electrodes
- 10 Ambu crocodile clips
- Protective hard carry case

\*On-going license costs apply

\*\*Terms and conditions apply

\*\*\*Additional costs may apply for on-site visits



For more information about the MAC 600 ECG System, visit our website: [www.numed.co.uk](http://www.numed.co.uk)  
call: **0114 243 3896** or email: [sales@numed.co.uk](mailto:sales@numed.co.uk)



SUPPLY | INSTALLATION | TRAINING | SUPPORT