

SpiroUSB

PC Based Spirometer

SpiroUSB is a spirometer that plugs directly into the USB port of your computer. It utilises the SPCS software allowing detailed spirometry testing to be performed.

Features and Benefits

- Supplied with intuitive SPCS software
- No batteries required – spirometer is powered directly via PC
- Highly portable spirometry system when used with a laptop
- Fully integrated with EMIS (additional 'Link Licence' required)
- Fully integrated with SystemOne
- Can upgrade existing Spida5 and SpidaXpert databases to SPCS
- 'Real time' Flow/Volume and Volume/Time graphs
- Test quality assessment to ensure accurate results
- Choice of four animated incentives including 'The Bubblegum Kid' and 'Wolf'
- Automatic predicted values calculation with ethnic correction factors
- 'Lung Age' estimation
- Reversibility testing with % change calculation
- Trend graphs can be created to show change in lung function over time
- Configurable printout format with Flow/Volume and Volume/Time curves
- NICE 2010 interpretation
- Calibration verification feature



PC Not Included

Leasing options available. Please enquire

2 Year Warranty

MediServe Support

- 90 Days Included
- 1, 2 or 3 Year Contracts Available

Optional Extras

- 3 Litre Accuracy Check Syringe

Consumables

- One Way Safety Mouthpieces
- PeraSafe Sterilising Agent

Fully Integrates with



Installation & Training Available