

NON-CONTACT INFRARED THERMOMETER

JPD-FR203



Professional
Accuracy



1 second
Reading



Hygienic
Non-Touch



Smart
Fever Alarm



Silent Mode



Child Mode

Features:

- Reliable and Accurate - Clinically tested and proven to be reliable and accurate. The most recommended infrared thermometer by doctors and pediatricians. Unique child temperature mode for more accurate measurement.
- Multifunctional- Measurement of forehead, objects & ambient temperature.
- Readings Retrieval - Can store and recall up to 20 readings.
- Sound control - Silent operation especially for measurements at night.
- Safe and Hygienic - No-touch option helps minimize spreading of germs. Completely safe for use on children and adults.

Specifications :

Display: HD LCD Display

Measuring Range & Accuracy:

Forehead mode: 35.0°C–42.2°C (95.0°F–108.0°F) ; $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$)

Object mode: 0.0°C–100.0°C (32.0°F–212.0°F); $\pm 1.0^{\circ}\text{C}$ ($\pm 2.0^{\circ}\text{F}$)

Battery: AAAx2

Memory Storage: 20 sets of measurement data

Temperature Units: °C/°F adjustable

Automatic Shutdown Function: 10s \pm 1s

Display Resolution : 0.1°C / 0.1°F



JUMPER

SHENZHEN JUMPER MEDICAL EQUIPMENT CO., LTD.

proactmedical.co.uk

Tel: +44 (0)1536 461981

E-mail: orders@proactmedical.co.uk

PROACT Medical Ltd

9-13 Oakley Hay Lodge
Great Folds Road
Corby
NN18 9AS