

# Patient Monitoring

## Vital signs monitors product comparison matrix



Portrait™ VSM Monitor



CARESCAPE™ VC150 Monitor



CARESCAPE V100 Monitor

Technical features			
Screen size	10.1 in (diagonal)	8.4 in (diagonal)	7.6 in (diagonal)
Dimensions	(H x W x D) 27.5 x 26.5 x 17.5 cm (10.8 x 10.4 x 6.8 in)	(H x W x D) 24.7 x 24.2 x 13.6 cm (9.75 x 9.5 x 5.3 in)	(H x W x D) 19.5 x 21.9 x 13.5 cm (7.7 x 8.6 x 5.3 in)
Waveform display	Yes	Yes	N/A
Weight	8.4 lb w/ battery	6.2 lbs w/ battery	5.4 lbs w/ battery
Touch screen type	Capacitive	Resistive	N/A
Printing capability	Yes, thermal dot ray	Yes, thermal dot ray	Yes, thermal dot ray
Wireless (WLAN)	Yes	Yes	No
HIS connectivity	Yes, through HL7 interface engine, or direct to certain IHE-compliant systems	Yes, through HL7 interface engine, or direct to certain IHE-compliant systems	No
RFID	Yes	No	No
Battery	Lithium Ion, 5.5 hrs operating time	Lithium-ion, 10.8V and 11.1V; up to 10.5 hrs operating time	Sealed lead acid, 6V up to 11.5 hrs operating time
Clinical features			
Continuous monitoring capability	Yes	Yes	Yes
NIBP technology	DINAMAP™	DINAMAP	DINAMAP
Pulse oximetry options	GE HealthCare's TruSignal SpO <sub>2</sub> , perfusion index Nellcor SpO <sub>2</sub> , SatSeconds Masimo SET SpO <sub>2</sub> , perfusion index	GE HealthCare's TruSignal SpO <sub>2</sub> , perfusion index Nellcor SpO <sub>2</sub> , perfusion index; SatSeconds, Saturation pattern detection (SPD) Masimo rainbow SET SpO <sub>2</sub> , perfusion index, SpHB, SpMet, SpCO, RRa, SpOC	GE HealthCare's TruSignal Nellcor OxiMax Masimo SET



	Portrait™ VSM Monitor	CARESCAPE™ VC150 Monitor	CARESCAPE V100 Monitor
<b>Respiration rate options</b>	Manual input. Timer feature supports accurate measurement For patients using the Portrait Mobile wireless monitoring solution, scan the QR code on the Portrait Mobile Patient Monitor, admit the patient to the Portrait VSM monitor, and the patient's RR (and PR and SpO <sub>2</sub> ) appear on the screen	RRa only available with the Masimo Rainbow SET technology with acoustic sensors Manual input as an observation item	N/A
<b>Temperature options</b>	Welch Allyn SureTemp, Exergen TemporalScanner; HeTaiDa non-contact infrared body thermometer	Welch Allyn SureTemp, Exergen TemporalScanner	Exergen TemporalScanner
<b>Round Advisor™ feature</b>	Yes	No	No
<b>Snapshots capability</b>	Yes, stores up to 10,000 entries	Yes, stores up to 3000 entries	Yes, stores up to 40
<b>Portrait Mobile monitoring solution compatibility</b>	Yes	No	No
<b>EWS capability, protocols</b>	Yes, 7 (NEWS+MEWS+5 customized)	Yes, 3	No
<b>Alarm event history</b>	Yes	No	No
<b>Automatically send the patient data to EMR in continuous monitoring mode</b>	Yes	No	No
<b>Lock screen for cleaning and disinfections</b>	Yes	No	No

Please refer to the product documentation for detailed information regarding configurations, specifications and capabilities.

Not all products or features are available in all markets. Full product technical specification is available upon request. Contact a GE HealthCare Representative for more information.

Data subject to change.

© 2024 GE HealthCare.

CARESCAPE, DINAMAP and Portrait are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

JB27767XX 01/2024