

# Compact trolley

## Instrument type:

Compact trolley designed to provide mobility for GE MAC 800 resting ECG system.

## Dimensions (Without resting ECG system):

Height: 37 in (94 cm)

Height with acquisition module arm and holder: 54 in (137 cm)

Width: 20 in (52 cm)

Width with barcode shelf: 22 in (56 cm)

Depth: 22 in (57 cm)

Depth with printer shelf: 25 in (65 cm)

Basic trolley weight: 31 lbs (14 kg)

Acquisition module arm and holder + Barcode shelf weight: 7 lbs (3 kg)

Printer shelf weight: 15 lbs (7 kg)

## Supported ECG System:

MAC 800



### Optional:

Printer Shelf

Simple Arm and Barcode kit (simple arm + barcode holder + card shelf)

Simple arm is compatible with value cable and multilink cable



Holder with multilink cable



Holder with Value leadwires cable

Barcode holder supports both new and old barcode scanner



Compact trolley is produced by SUZHOU TIANYE  
MECHANICAL INDUSTRY CO.,LTD for GE Healthcare  
MAC800 resting ECG system. SUZHOU TIANYE  
MECHANICAL INDUSTRY CO.,LTD owns the compact  
trolley legally.